

The IOT Based Jacket

Eng. Elie El Arou
Eng. Atieh Atieh



Outline



Introduction

What is the problem



Our Solution!

Requirements
System Design
Functional Decomposition



Components

The hardware and
specification



Hardware Implementation

Circuit and flow
chart



Software Implementation

Flow Charts, Code Snips
And application



Project Analytics

Cost, Gantt Chart, Risk
Assessment, and Conclusion



INTRODUCTION



STATISTICS (1)



ONE IN FOURTH

one in fourth of the elderly population report falling each year



300,000 HOSPITALIZED

There is 300,000 elderly individuals hospitalized annually due to broken hips



\$ 6,000,000,000 COSTS

medical costs for broken hips in the U.S. is six billion dollars per year

STATISTICS (2)

33%

33% of elderly individuals do not survive beyond one year after experiencing a hip fracture



31.8%

31.8% of elderly mortality in the US is caused by heart disease

WE HAVE THE SOLUTION!



PROJECT REQUIREMENTS

ELECTROCARDIOGRAM

Monitor live ECG

O2 LEVEL & HEART RATE

Monitor heart rate and oxygen level with 85 to 90% accuracy

TEMPERATURE

Track body and room temperature with 98% accuracy



FALL DETECTION

Track any fall and deploy airbag


WARMING MODULE

Give warmth through heated sheets


INTERFACE

Displayed and send notification on a web and PWA application

FUNCTIONAL DECOMPOSITION LEVEL 0



FUNCTIONAL DECOMPOSITION LEVEL 1



HARDWARE COMPONENTS

ELECTROCARDIOGRAM

AD-8232



MICROCONTROLLERS

ESP-8266 & ESP-WROOM-32



BLOOD OXYGEN LEVEL

MAX-30100



FALL DETECTION

MPU-6050



BODY TEMPERATURE

DS18B20





Heating pads

Heating pads


ESP-8266 Specifications

- Price: 3.5\$
- Microcontroller: Tensilica 32-bit RISC CPU Xtensa LX106
- Operating Voltage: 3.3V
- Digital I/O Pins (DIO): 16
- Analog Input Pins (ADC): 1
- UARTs: 1
- SPIs: 1
- I2Cs: 1
- Flash Memory: 4 MB
- SRAM: 64 KB
- Clock Speed: 80 MHz
- PCB Antenna



ESP-WROOM-32


Specifications



- Microprocessor: Tensilica Xtensa LX6
- Maximum Operating Frequency: 240MHz
- Operating Voltage: 3.3V
- Analog Input Pins: 12-bit, 18 Channel
- DAC Pins: 8-bit, 2 Channel
- Digital I/O Pins: 39
- DC Current on I/O Pins: 40 mA
- DC Current on 3.3V Pin: 50 mA
- SRAM: 520 KB
- Communication: SPI(4), I2C(2), I2S(2), CAN, UART(3)
- Wi-Fi: 802.11 b/g/n
- Bluetooth: V4.2 – Supports BLE and Classic Bluetooth

AD-8232

Specifications





- Price: 11\$
- Fully integrated single-lead ECG front end
- Supply current: 170 μ A
- Common-mode rejection ratio: 80 dB
- Two or three electrode configurations
- Accepts up to ± 300 mV of half cell potential
- Fast restore feature improves filter settling
- Leads off detection: ac or dc options
- Integrated right leg drive (RLD) amplifier
- Single-supply operation: 2.0 V to 3.5 V

AD-8232 Specifications


PINS:

- SDN
- LO+
- LO-
- OUTPUT
- 3.3V
- GND
- right arm (RA)
- left arm (LA)
- and right leg (RL)



MAX30100

Specifications





- Two light-emitting diodes
- Two photo detectors, optimized optics and low noise simulation.
- Communication modes: standard IIC protocol
- Power supply: 3.3V to 5.5V
- Current draw: ~600 μ A (during measurements)
- Current draw: ~0.7 μ A (during standby mode)
- Red LED Wavelength: 660nm
- IR LED Wavelength: 880nm
- Temperature Range: -40°C to +85°C
- Temperature Accuracy: $\pm 1^\circ\text{C}$


MAX-30100 Specifications

PINS:

- VCC
- GND
- SCL
- SDA
- INT
- IRD
- RD



MPU-6050 Specifications





- Power supply :3-5v
- Communication modes: standard IIC communications protocol
- DEGREES OF FREEDOM : 6 x
- Chip built-in 16bit AD converter, 16-bit data output
- Acceleration range: $\pm 2 \pm 4 \pm 8 \pm 16g$
- Pin spacing 2.54mm

MPU-6050 Specifications

PINS:

- VCC
- GND
- SCL
- SDA
- XCL
- XDA
- ADD
- INT



DS18B20 Specifications




- 1-Wire interface requires only one port pin for communication
- Requires no external components
- Power supply range is 3.0V to 5.5V
- Zero standby power required
- Measures temperatures from -55°C to $+125^{\circ}\text{C}$ (-67°F to $+257^{\circ}\text{F}$)
- $\pm 0.5^{\circ}\text{C}$ accuracy from -10°C to $+85^{\circ}\text{C}$

DS18B20 Specifications

PINS:

- VCC
- GND
- DATA




AD-8232

Specifications




- Material: Silicon Rubber
- Dimensions: 5x5cm
- Voltage input: 12VDC
- Power consumption: 10W

CIRCUIT IMPLEMENTATION



ARDUINO FLOW CHART



ARDUINO SETUP

Wi-Fi connection:

```
// Connect to WiFi
WiFi.begin(ssid, password);
while (WiFi.status() != WL_CONNECTED) {
  delay(2000);
  Serial.println("Connecting to WiFi...");
}
Serial.println("Connected to WiFi");
Serial.println("WiFi connected");
Serial.println("IP address: ");
Serial.println(WiFi.localIP()); //You can get
server.listen(80);
Serial.print("Is server live? ");
Serial.println(server.available());
```

MPU-Gyroscope:

```
if (!mpu.begin())
{
  Serial.println("Failed to find MPU6050 chip");
  while (1)
  {
    delay(10);
  }
}
mpu.setAccelerometerRange(MPU6050_RANGE_8_G);
```


ARDUINO SETUP (2)

MAX-30100 oxygen level:

AD-8232 ECG

```
Serial.print("Initializing pulse oximeter..");  
if (!pox.begin()) {  
  Serial.println("FAILED");  
  for (;;);  
} else {  
  Serial.println("SUCCESS");  
}  
// Configure sensor to use 7.6mA for LED drive  
pox.setIRLedCurrent(MAX30100_LED_CURR_7_6MA);  
// Register a callback routine  
pox.setOnBeatDetectedCallback(onBeatDetected);  
Serial.println("");
```

```
pinMode(4, INPUT);  
pinMode(2, INPUT);
```

OXYGEN LEVEL MAX-30100


```

pox.update();
if (millis() - tsLastReport > REPORTING_PERIOD_MS) {
    Serial.print("bpm / SpO2:");
    Serial.print(pox.getSpO2());
    Serial.println("%");

    // Prepare the JSON data to send
    String jsonData = String(String(pox.getSpO2()));
    if (client.available()) {
        client.send(jsonData);
        maxim.resetFifo();
    }

    tsLastReport = millis();
}

```



FALL DETECTION MPU-6050

```

if (millis() - lastMPUReadTime >= 500) {
    sensors_event_t a, g, temp;
    mpu.getEvent(&a, &g, &temp);
    roomT = "temperature2=" + String(temp.temperature);
    float accelMagnitude = (sqrt(a.acceleration.x * a.acceleration.x + a.acceleration.y * a.acceleration.y + a.acceleration.z * a.acceleration.z)) - 10;


    if (accelMagnitude > 6)
    {
        // Potential fall detected
        Serial.println("Fall detected");
        FallDetected();
    }
}

```

```

void FallDetected()
{
    HTTPClient https;
    std::unique_ptr<BearSSL::WiFiClientSecure> client(new BearSSL::WiFiClientSecure);
    client->setInsecure();
    Serial.print("[HTTPS] begin...\n");
    if (https.begin(*client, serverUrl)) {
        Serial.print("[HTTPS] POST...\n");
        https.addHeader("Content-Type", "application/json");
        // Create a JSON document and populate it with your data
        DynamicJsonDocument jsonDoc(256);
        jsonDoc["fallstatus"] = 1;
        // Serialize the JSON document to a string
        String jsonPayload;
        serializeJson(jsonDoc, jsonPayload);
        int httpCode = https.POST(jsonPayload);
        if (httpCode > 0) {
            Serial.printf("[HTTPS] POST... code: %d\n", httpCode);
            if (httpCode == HTTP_CODE_OK || httpCode == HTTP_CODE_MOVED_PERMANENTLY) {
                String payload = https.getString();
                Serial.println(payload);
            }
        } else {
            Serial.printf("[HTTPS] POST... failed, error: %s\n", https.errorToString(httpCode).c_str());
        }
        https.end();
    } else {
        Serial.printf("[HTTPS] Unable to connect\n");
    }
}

```



ECG AND HEART RATE AD-8232

```

unsigned long lastPeakTime = 0;
unsigned long currentPeakTime = 0;
int threshold = 400;
bool peakDetected = false;
float bpm = 0;


void loop() {
  delay(20);
  unsigned long currentTime = millis();
  Serial.println(analogRead(34));
  int ecgValue = analogRead(34);
  if (ecgValue > threshold && !peakDetected) {
    peakDetected = true;
    currentPeakTime = millis();
    if (lastPeakTime != 0) {
      float time = currentPeakTime - lastPeakTime;
      // Calculate heart rate in beats per minute (BPM)
      if( time > 400 ){
        bpm = 60000.0 / time ;
        Serial.println(bpm);
      }
    }
    lastPeakTime = currentPeakTime;
  }
  if (ecgValue < threshold) {
    peakDetected = false;
  }
}

```

```

String JsonData = "{\"timestamp\":" + String(currentTime) + ", \"ecgVal\":" + String(ecgValue) + ", \"heartRate\":" + String(bpm)}";
if (WiFi.status() == WL_CONNECTED)
{
  server.poll();
  if (!client.available()) {
    client = server.accept();
    if (client.available()) {
      Serial.println("Client connected");
    }
  }
  if (client.available()){
    client.send(JsonData);
  }
}

```




TEMPERATURE SENSOR DS18B20

```


unsigned long currentMillis = millis();
if (currentMillis - lastTemperatureUpdateTime >= temperatureUpdateInterval)
{
    // It's time to update temperature
    sensors.requestTemperatures();
    float temperature = sensors.getTempCByIndex(0);
    sendTemperatureToServer(temperature);
    lastTemperatureUpdateTime = currentMillis; // Reset the timer
}

void sendTemperatureToServer(float temperature)
{
    String data = "temperature=" + String(temperature) + "&" + roomT;
    std::unique_ptr<BearSSL::WiFiClientSecure> client(new BearSSL::WiFiClientSecure);
    client->setInsecure();
    HTTPClient https;
    Serial.print("[HTTPS] begin...\n");
    if (https.begin(*client, "https://JackBack.onrender.com/api/update")) {
        Serial.print("[HTTPS] POST...\n");
        https.addHeader("Content-Type", "application/x-www-form-urlencoded");
        int httpCode = https.POST(data);
        if (httpCode > 0) {
            Serial.printf("[HTTPS] POST... code: %d\n", httpCode);
            if (httpCode == HTTP_CODE_OK || httpCode == HTTP_CODE_MOVED_PERMANENTLY)
            {
                String payload = https.getString();
                Serial.println(payload);
            }
        }
    }
}

```



HARDWARE TESTING



FALL DETECTION MPU-6050

Natural Movements

-9.33	0.51	0.88	0.03	0.27	0.28	1.36	0.71	1.46	1.01	1.15
0.76	0.78	0.84	0.72	0.91	0.73	0.81	0.75	0.29	0.82	0.54
0.39	0.69	1.16	0.63	0.63	1.21	0.53	1.61	0.52	1.11	0.64
1.02	0.69	-0.21	1.47	1.29	0.33	-0.16	1.26	1.77	0.92	0.52
1.23	0.35	0.70	0.94	1.02	0.43	0.89	1.18	0.70	1.46	0.68
1.21	1.30	1.12	1.48	-0.17	1.38	-0.02	2.17	-0.11	1.31	0.47
0.58	0.06	0.34	2.36	1.79	0.69	1.60	1.17	1.11	-0.86	0.48
-0.45	1.22	0.93	0.25							

FALL DETECTION MPU-6050

Fall Detection with threshold [-2,+2]

0.80	0.86	2.04	Fall detected	
0.96	0.05	1.27	1.33	2.11 Fall detected
0.72	0.13	23.85	Fall detected	
2.39	Fall detected			
0.33	0.45	0.71	2.10	Fall detected
0.62	0.72	0.42	0.83	0.56 1.14 .144 0.53 1.76 1.93
0.47	2.46	Fall detected		
0.76	1.47	0.88	0.80	-2.50 Fall detected
0.99	1.28	1.65	-5.27	Fall detected
0.69	-0.62	1.60	-0.13	-2.18 Fall detected
0.10	0.16	0.41	0.42	-3.55 Fall detected
1.11	1.10	1.11	1.01	1.08 1.01 1.15 1.09 0.99 1.13
1.07	9.84	Fall detected		
0.78	0.96	1.23		

FALL DETECTION MPU-6050

```

0.39  0.43
0.37  0.41  0.34  0.40  0.41  0.40  0.42  0.43  0.39  0.41  -0.27
-0.47 0.43  0.90  -0.59  30.74 Fall detected
0.11  -0.40  -4.33 Fall detected
-2.07 -1.40 16.54 Fall detected
-0.48  4.90 Fall detected
-0.48 0.04  -1.62 -0.62 -0.50  4.45 Fall detected
-3.56 -1.68 -1.17 0.40  0.29

```

Fall Detection with threshold [-4,+4]

Fall Detection with threshold [-6,+6]

```

                                1.05 0.88 1.21 1.30 0.81 0.07 11.73 Fall detected
                                0.26 1.93 2.26 -0.74 0.84 1.07 -6.54 Fall detected
                                1.18 0.79 0.53 0.91 0.50 10.86 Fall detected
                                1.30 2.13 1.87 1.14 1.76 2.60 1.52 2.56 2.07 -0.22 1.92
                                3.99 0.96 2.31 2.71 2.81 3.85 1.12 0.85 2.77 1.30 1.36
                                3.77 0.70 2.08 1.15 7.59 Fall detected
                                1.17 0.98 1.06 1.70 2.16 0.59 0.61 0.44 9.90 Fall detected 0.06
                                4.93 -1.20 -0.64 19.29 Fall detected
                                                -6.84 Fall detected
                                                0.64

```

AIRBAG



1

Fall Detection Delay

500 milliseconds

2

Airbag Delay

50 milliseconds

3

Total Delay

550 milliseconds

4

Equation

$$time = \sqrt{\frac{2h}{g}} = 580 \text{ milliseconds}$$

HEATING PADS



1

Standard temperature

According to Cardinal Health, the appropriate and safe temperature range is 40-45°C.

2

At 9 volts



Temperature is 35.63°C

3


At 12 volts

Temperature is 43.94°C


ECG AND HEART RATE AD-8232




OXYGEN LEVEL MAX-30100




1% error







APPLICATION CHART



IN DEPTH APPLICATION CHART





LIVE MONITOR

Temperature hide ^

Body Temperature: 37 °C

Room Temperature: 20 °C

HeartRate show v

OxygenLevel show v

Google Chrome




Room Temperature Hazard!

Room temperature is below 21 it is 20

all-in-one-jacket.web.app

TEMPERATURE

Arduino Side



Node.js Side

```
void sendTemperatureToServer(float temperature)
{
    String data = "temperature=" + String(temperature) + "&" + roomT;
    std::unique_ptr<BearSSL::WiFiClientSecure> client(new BearSSL::WiFiClientSecure);
    client->setInsecure();
    HTTPClient https;
    Serial.print("[HTTPS] begin...\n");
    if (https.begin(*client, "https://JackBack.onrender.com/api/update")) {
        Serial.print("[HTTPS] POST...\n");
        https.addHeader("Content-Type", "application/x-www-form-urlencoded");
        int httpCode = https.POST(data);
        if (httpCode > 0) {
            Serial.printf("[HTTPS] POST... code: %d\n", httpCode);
            if (httpCode == HTTP_CODE_OK || httpCode == HTTP_CODE_MOVED_PERMANENTLY)
```

```
app.post('/api/update', async (req, res) => {
    console.log(req.body);
    const bodytemperature = req.body.temperature;
    const roomtemperature = req.body.temperature2;

    io.sockets.emit('TempUpdate', { bodytemperature, roomtemperature });
    console.log(parseInt(roomtemperature));

    if (parseInt(roomtemperature) < 21) {
        console.log('reach');
        const TempnotificationPayload = {
            notification: {
                title: 'Room Temperature Hazard!',
                body: `Room temperature is below 21 it is ${roomtemperature}`
            }
        };

        try {
            const subscribers = await fetchSubscribersFromDatabase();
            await Promise.all(subscribers.map(sub => webpush.sendNotification(sub, JSON.stringify(TempnotificationPayload))));
```

TEMPERATURE

Temperature Service



Subscribe

Home-Page Component

```
import { Observable } from 'rxjs';

@Injectable({
  providedIn: 'root'
})
export class TemperatureService {

  constructor(private http: HttpClient, private socket: Socket) {}

  getTemperatureSocket(): Observable<any> {
    return this.socket.fromEvent('TempUpdate');
  }

  getLastTemperature(): any {
    const storedTemperature = sessionStorage.getItem('temperature');
    return storedTemperature ? JSON.parse(storedTemperature) : '';
  }

  setLastTemperature(data: any): void {
    sessionStorage.setItem('temperature', JSON.stringify(data));
  }
}
```

```
this.Tempsocket = this.TempService.getLastTemperature() ;

this.TempService.getTemperatureSocket().subscribe((message: any) => {
  this.TempService.setLastTemperature(message);
  this.Tempsocket = message
  console.log(this.Tempsocket)
  console.log(message)
});
```

Html Template Data Binding

```
> Body Temperature: {{Tempsocket?.bodytemperature}} °C </div> <br>
```



FALL DATA + NOTIFICATION

```
app.post('/api/FallDetected', async (req, res) => {
  const fallStatus = req.body.fallstatus;
  const currentDate = new Date();
  const formattedTime = currentDate.toLocaleString('en-US', {
    timeZone: 'EET'
  });
  const formattedDate = currentDate.toDateString();
  const IsoData = formattedTime.split(',');
  const Isod = currentDate.toISOString().split('T')[0];
  const IsoTime = IsoData[1]
  console.log(formattedDate+ " " + formattedTime)

  const FallEvent = {
    IsoDate: Isod,
    Date: formattedDate,
    Time: IsoTime
  }
  try{ const response = await db.collection("Falls").add(FallEvent);
}
catch(err){
  console.log(err)
}
```

```
try {
  const subscribers = await fetchSubscribersFromDatabase();
  await Promise.all(subscribers.map(sub => webpush.sendNotification(sub, JSON.stringify(notificationPayload))))
  res.status(200).json({ message: 'Notifications sent successfully.' });
} catch (err) {
  console.error("Error sending notifications:", err);
  res.sendStatus(500);
}
```

```
requestNotificationPermission() {
  Notification.requestPermission().then(permission => {
    if (permission === 'granted') {
      this.subscribeToNotifications();
    }
  });
}
```



FALL DATA QUERY

Start date: 11/10/2023

End date: 11/11/2023

NOV 2023

Fall Date	Fall Time
Sat Nov 11 2023	6:51:21 PM
Sat Nov 11 2023	1:54:21 PM
Sat Nov 11 2023	1:54:27 PM
Sat Nov 11 2023	6:05:27 PM

```
app.get('/api/ReadFall/:StartDate/:EndDate', async (req, res) => {
  const start = req.params.StartDate;
  console.log(start);
  const end = req.params.EndDate;
  console.log(end);
  const query = db.collection('Falls')
    .where('IsoDate', '>=', start)
    .where('IsoDate', '<=', end);

  try {
    const getQuery = await query.get();
    const falls = [];

    getQuery.forEach((doc) => {
      falls.push(doc.data());
    });

    res.json(falls);
  } catch (err) {
    console.error(err);
    res.send(err);
  }
});
```



ECG PLOTTING

```
ws2.on('open', () => {
  console.log('Connected to WebSocket server');
  ws2.send('Hello from the client!');
});

var dataArray = [];

ws2.on('message', (Data) => {
  const decodedString = Data.toString('utf-8');
  // console.log('Received message:', decodedString);
  dataArray.push(decodedString);
  if(dataArray.length >= 50)
  {
    dataArray = dataArray.slice(dataArray.length-50)
  }
});

setInterval(() => {
  io.sockets.emit('ECG', dataArray.slice()); // Send a copy to prevent modification issues
}, 2500);
```




ECG SAVING

```
app.post('/api/SaveECG', async (req, res) => {
  const { name, points } = req.body;

  console.log('Received data:', { name, points });

  try {
    const existingDoc = await db.collection("ECG").doc(name).get();





    if (existingDoc.exists) {
      console.log('yes')
      res.json({ success: false, error: 'Document with the same name already exists' });
    } else {
      await db.collection("ECG").doc(name).set({ points });
      res.status(200).json({ success: true, chartId: name });
    }
  } catch (err) {
    console.error(err);
    res.status(500).json({ success: false, error: 'Internal Server Error' });
  }
});
```



ECG PLOTTING & SAVING

```
app.get('/api/GetAllSavedECGNames', async (req, res) => {
  console.log('namesReached')
  try {
    const usersRef = db.collection("ECG");
    const response = await usersRef.get();

    const docNames = response.docs.map(doc => doc.id);
    res.send(docNames);
  } catch (err) {
    res.status(500).send(err);
  }
});
```



```
app.get('/api/GetSavedECG', async (req, res) => {
  try {
    const chartName = req.query.name;
    console.log(chartName)

    const userRef = db.collection("ECG").doc(chartName);
    const doc = await userRef.get();
    if (doc.exists) {
      res.send(doc.data().points);
    } else {
      res.send({ error: "Document not found" });
    }
  } catch (err) {
    res.status(500).send(err);
  }
});
```

ANGULAR PWA




- Modern web technology to make websites feel more like the apps.
- Use functionalities that are only for application on a website like notification.
- Market of PWA will reach up to 10.77 billion dollars by 2027
- Google reports a significant increase of 270% in desktop installations of PWAs from the beginning of 2021.



TO INSTALL THE APP SCAN THE QR CODE




RISK ASSESSMENT USING FMEA



Types of Rating	Heater	ECG	Oxygen	Airbag	Software
Severity Ratings	9	6	3	10	7.5
Occurrence Ratings	5	1	1	6.5	5.5
Detection Ratings	8	8	8	8	6
RPN	360	48	48	520	247.5

COMPONENTS	COST IN DOLLARS
ESP-8266	\$3.50
ESP-WROOM-32	\$7.77
AD-8232 + ELECTRODES	\$11.00 + \$6.00
MAX-30100	\$7.00
MPU-6050	\$4.00
DS18B20	\$2.00
HEATING PADS	\$6.00 * 2
JACKET	\$15.00
POWER BANK	\$20.00
ADDITIONAL MATERIAL (BREADBOARD, WIRES ETC..)	\$20.00
TOTAL COST	\$108.27

WORK BREAKDOWN



Activity	JUL	AUG	SEP	OCT	NOV
GET ALL COMPONENTS NEEDED	ELIE & ATIEH	ELIE & ATIEH			
PROGRAM THE HEART RATE			ELIE		
PROGRAM THE OXYGEN SENSOR			ATIEH		
PROGRAM THE TEMPERATURE SENSOR			ELIE		
PROGRAM FALL SENSOR			ATIEH		
DESIGN JACKET				ELIE & ATIEH	
ADD ALL COMPONENTS				ELIE & ATIEH	
CREATE THE FRONT-END				ELIE & ATIEH	
CREATE THE BACK-END				ELIE & ATIEH	
DEPLOY ON SERVER					ELIE
TEST ALL COMPONENTS					ELIE & ATIEH

THANK YOU!



KEEP YOUR ELDERLY SAFE